

Agilent BioTek 406 FX Microplate Washer Dispenser **User Manual**



Notices

Document Identification

2121000N

Revision C. TBD 2024

Copyright

© Agilent Technologies, Inc. 2021-2023

No part of this manual may be reproduced in any form or by any means (including electronic storage and retrieval or translation into a foreign language) without prior agreement and written consent from Agilent Technologies, Inc. as governed by United States and international copyright laws.

Agilent Technologies, Inc. 5301 Stevens Creek Blvd. Santa Clara, CA 95051 USA

Software Revision

This document is valid for the 3.0 revision or higher of the Agilent Liquid Handling Control (LHC) and compatible programs until superseded.

Instrument Manufacturing



Manufactured by Agilent Technologies, Inc. 5301 Stevens Creek Blvd. Santa Clara, CA 95051



Warranty

The material contained in this document is provided "as is," and is subject to being changed, without notice, in future editions. Further, to the maximum extent permitted by applicable law, Agilent disclaims all warranties, either express or implied, with regard to this manual and any information contained herein, including but not limited to the implied warranties of merchantability and fitness for a particular purpose. Agilent shall not be liable for errors or for incidental or consequential damages in connection with the furnishing, use, or performance of this document or of any information contained herein. Should Agilent and the user have a separate written agreement with warranty terms covering the material in this document that conflict with these terms, the warranty terms in the separate agreement shall control.

Technology Licenses

The hardware and/or software described in this document are furnished under a license and may be used or copied only in accordance with the terms of such license.

Restricted Rights Legend

U.S. Government Restricted Rights, Software and technical data rights granted to the federal government include only those rights customarily provided to end user customers.

Agilent provides this customary commercial license in Software and technical data pursuant to FAR 12.211 (Technical Data) and 12.212 (Computer Software) and, for the Department of Defense, DFARS 252.227-7015 (Technical Data - Commercial Items) and DFARS 227.7202-3 (Rights in Commercial Computer Software Documentation).

Safety Notices

CAUTION

A **CAUTION** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in damage to the product or loss of important data. Do not proceed beyond a **CAUTION** notice until the indicated conditions are fully understood and met.

WARNING

A WARNING notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in personal injury or death. Do not proceed beyond a WARNING notice until the indicated conditions are fully understood and met.

Contents

No	otices	
1	General Information	1
	Contact Information	2
	Worldwide Sales and Support	2
	Technical Support and Service	2
	Customer Care	2
	Intended Use Statement	2
	Quality Control	2
	Safety Notices	2
	Warnings and Precautions	
	Electrical Hazards	3
	Chemical/Environmental	4
	Components	5
	Intended Product Use	6
	Safety Symbols	6
	Conformance to Standards and Directives	12
	Directives	12
	Standards	12
	EMC Information and Technical Description	13
	Emissions	
	Immunity	13
	Disposal	13
2	Introduction	14
	Introducing the 406 FX Microplate Washer Dispenser	15
	Package Contents	16
	Accessory Pack	16
	Optional Accessories	16
	General Instrument Accessories	16
	Washer Accessories	17
	Washer Manifolds	18
	Magnetic Bead Assay Accessories	18
	Miscellaneous	18
	Peri-pump Optional Accessories	19
	Peri-pump Modules	19
	Cassettes and Accessories	19
	Tips and Cassette Tubing Extension Kits	20
	General accessories	20

	Syringe Dispenser Optional Accessories	21
	Manifolds	21
	Bottles and Tubing	21
	General accessories	21
	Integration Kits	22
	Features of the 406 FX	22
	Washer	22
	Peri-pump Dispenser	23
	Syringe Dispensers	23
	Liquid Handling Control (LHC) Software	23
	BioStack Compatibility	24
	406 FX Dispenser Comparison	25
3	Installation	27
	Important Information	
	Unpack the Instrument	
	Remove the Shipping Hardware	
	Supplies	
	Setting Up the 406 FX	
	Connect the Vacuum Pump, Tubes, and Bottles	
	Waste System	31
	Install the Vacuum Line Filter	33
	Final Check	33
	Attach the Mist Shield	33
	Install the Microplate Carrier	34
	Install the Left-Side Carrier Kit (optional)	35
	Set Up the Peri-pump Dispensers	36
	Install the Secondary Peri-Pump	36
	Dispense Cassette Diagram	39
	Install the Dispense Cassette	40
	Change the cassette type setting	42
	Prime Trough Inserts	42
	Install the Syringe Dispenser Component	44
	Install Inline Filter for 32-Tube Dispensers	45
	Install Tubing and Manifolds for Syringe Dispenser	46
	Check the Syringe Manifold Installation	50
	Syringe Dispenser Check Valves	50
	Install Software/Connect to Computer	50
	Connect to Host Computer	51
	Connect Power Cable	51

	Define Instrument Settings	51
	LHC Users Only	51
	Standard Vacuum Pump Users	52
	Define Startup Preferences (LHC users)	52
	Verify Performance	52
	Verify the Peri-Pump Dispenser	53
	Verify the Syringe Dispensers	53
	Verify the Washer Component	54
	Repackaging and Shipping Instructions	55
	Repacking the Syringe Dispenser	61
4	Operation	63
	Basic Operation	
	Optimize Performance	
	Keep the devices clean and the tubing wet	
	Prime the tubing to remove air bubbles	
	Best Practices for all 406 FX Devices	
	Before using the Peri-pump	65
	Before using the Syringe Dispenser	
	Optimize protocols to improve evacuation	66
	Optimize protocols to specify the desired residual volume	66
	RUN: Running Predefined Protocols	
	Creating Protocols: Washing, Aspirating and Dispensing Fluid	67
	Touch Screen Control	67
	Liquid Handling Control (LHC) Software	67
	Customize the Predefined Protocols	68
	Predefined Protocols	69
	Predefined Sample Protocols	70
	Protocol Parameters Tables	72
	Syringe Dispenser Flow Rates	74
	Wash Step Parameters Table	75
	Operating with the BioStack	77
	Touch Screen Control:	77
	LHC Control:	78
	Touch Screen Basics	79
	Quick Prime	80
	Quick Prime	80
	Quick Wash	80
	Quick Wash	80
	Learn more:	80

Quick Clean	80
Quick Clean (Sonicate)	80
Add protocols to the Home screen	81
Remove protocols from the Home screen	81
Create, Edit, or Copy a Protocol	82
Copy a Protocol	83
Lock (password-protect) and unlock protocols	84
Rename a protocol	84
Change the Plate Type	85
Using LHC to Control the 406 FX Washer Dispenser	85
Predefined Protocols	85
Communications Port	85
406 FX Steps	86
LHC Protocols	86
Customize the Protocols	
LHC Protocols Explained	87
Prerequisite	87
Protocol Files	87
Managing Multiple Instruments	88
Validate versus Run	88
Target Instrument Settings	89
Get settings from:	89
How to define a Protocol (LHC only)	90
Customize the Predefined Protocols	91
How to define the Plate Map	91
Plate Types and Processing Patterns	92
Plate Types Table	93
Handling special plates	93
1536F - 1536-well Flanged Plates	93
Handling Special Plates	94
1536F - 1536-well Flanged Plates	94
About Wash Processing Patterns	95
Dispense Processing Patterns	97
Columns	97
Rows	97
Peri-pump Dispense Pattern	98
Dispense Pattern with Both Syringes	102
Washer Operation	102
About the 406 FX Wash Step	102

Wash Step Parameters Table	102
Shake/Soak Step Parameters	
Pre-Dispense	105
Washer Flow Rates	106
Recommendations for 1536-well hardware	107
How to define a Final Aspiration step	108
Add a secondary aspiration to a wash cycle	108
Cell Wash	
Low-flow Fluid Path	109
Define a Cell Wash Protocol	110
Cell Wash Strategies	
Best practices:	112
Biomagnetic Separation - Magnetic Bead Assays	113
Handling and Cleaning the Magnets	
Realign the BioStack with the Magnet Installed	114
Perform Magnetic Bead Assays	
Optimize Magnetic Bead Protocols	114
Flat Magnets	
Ring Magnets	117
How to determine the Magnet Adapter Height Offset	118
Magnet Height Offset	118
Bottom-Wash	119
Washer Flow Rates	120
Peri-pump Peristaltic Dispenser	121
Quick Dispense-Prime-Purge	121
Dispense Cassettes	122
Peri-pump: How it works	122
Recommended prime volumes for the Peri-pump	123
At the start of the day:	123
Release the tension on the dispense cassette	124
Peri-pump Dispense Step	125
Define the Dispense Step	125
Require a Specific Peri-pump Cassette	125
Shorten the Dispense Cassette Tubing	126
Standard Cassette Dead Volume	127
Peri-pump Advanced Settings	127
Settings for 0.5 µL Dispense	127
Dual Syringe Dispenser	128
Quick Dispansa	120

How to Prime the Syringe dispenser	128
Recommended prime volumes for the Syringe dispensers	130
Syringe Prime Step	130
Syringe Dispense Step	131
Change the Syringe Dispenser Manifold	132
Special Procedure for Magnetic Bead Assays	133
406 FX Instrument Utilities: General	135
Adjust Utility	137
Positional Limits	137
How to use the Adjust Utility- LHC	137
Use the Align Tool to optimize manifold tube positioning	138
AutoPrime	139
About AutoPrime	139
Specify the Interval and AutoPrime Parameters	
Turn on AutoPrime for the Washer	140
Turn on AutoPrime for the Peri-pump Dispenser	
Turn on AutoPrime for the Syringe Dispenser	
Delay After Dispense	
Standard Vacuum Pump Users	
Upload-Download Protocols (LHC Only)	
Operating with the BioStack	
Touch Screen Control:	
LHC Control:	
BioStack Alignment Utility	
Plate Lid Definition	
Transfer Plate Lid Definitions	
Use the BioStack	
Re-stack plates in the BioStack	
Nth Plate to perform special process	
Using Nth Plate to perform special functions	
About the Nth Plate Step	
Advanced Plate Sequencing with the BioStack	
Maintenance	
Preventative Maintenance Overview	
Washer	
Peri-pump	
Syringe Dispenser	
Recommended Maintenance Schedule	
Recommended Maintenance Schedule	160

5

Daily Maintenance	161
Overnight/Multi-Day Maintenance	162
AutoPrime for Overnight-Weekend Maintenance	162
Soak the manifold tubes in cleaning fluid	164
Prime Step	164
Create a Prime protocol	165
Removing Protein Residuals and Fungi Growth	166
Daily Practice with buffer and deionized water:	166
Weekly or As Needed, use NaOH and HCl to remove proteins:	166
Alternatively, use an Enzyme-Active Detergent:	167
Direct Drain Only: Clean the Waste Bottle Filter	167
Periodic Maintenance	
Clean the Bottles	
Clean the Plate Carrier	169
Clean the Exterior Surfaces and Mist shield	169
Autoclavable Components	169
Washer Maintenance	170
AutoClean the Washer	170
Create an AutoClean protocol	171
Touch screen	172
AutoClean the 32-Tube Dispense Manifold	172
Tools for cleaning dispense and aspirate tubes	173
Remove and clean the washer manifold	173
Clean the Buffer Bottle Filter	175
Peri-pump Dispenser Maintenance	175
Flush the Dispense Cassette	176
Tools and Supplies	176
Unclog the Dispense Tips	177
Record the Number of Plates Processed	178
Cassette Expected Lifetime	178
Syringe Dispenser Maintenance	178
Clean the Bottles and Tubes	179
Syringe Dispenser Maintenance Tasks	179
Clean the Manifolds	179
Clean the Dispense Tubes	180
Clean or Replace the Check Valves	181
Autoclave the Syringe Head	182
Decontamination	183
Tools and Supplies	103

	Decontaminate Exterior Surfaces	184
	Decontaminate Tubing and Manifold	184
	Alternate Decontamination Procedure for Tubing and Manifold	187
	Long Shutdown (Prepare for Storage or Shipment)	188
	Storing the Instrument	189
	Replace Components	189
	Replace Washer Manifold O-rings and Seals	189
	Replace Dispense Cassette Tips	191
	Calibrate the Backlash for Syringe Dispenser	192
	Plate-Type-Specific Syringe Manifolds	192
6	Qualification	194
	Qualification Overview	
	Washer	195
	Washer Qualification Tests	195
	Peri-pump and Syringe Dispensers	
	Qualification Schedule	
	Liquid Testing the 406 FX	196
	Which Tests to Perform?	196
	Washer	196
	Peri-pump Dispenser	197
	Syringe Dispenser	197
	Important Recommendations for All Liquid Tests	197
	System Self-Test, Verify Information	198
	Record Onboard Software Information	198
	Sensor Test: Fluid, Waste	198
	Test the Fluid Supply Sensor	198
	Test the Waste Level Sensor	199
	Ultrasonic Cleaner Test	199
	Washer Dispense Precision and Accuracy Test	199
	Predefined Test Protocols	200
	Predefined Protocols	200
	Washer Qualification Test Materials	201
	Evacuation Efficiency Test (96)	202
	Evacuation Diagnostic Test	
	Dispense Precision and Accuracy Test (96)	
	Buffer Switching Dispense Test	
	Evacuation Efficiency Test (192)	
	Evacuation Diagnostic Test (192)	
	Dispense Precision and Accuracy Test (192)	207

Evacuation Efficiency Test (128)	208
Dispense Precision and Accuracy Test Worksheet	
96-Tube Manifold	
Evacuation Efficiency Test Worksheet	210
96-Tube Manifold	210
Dispense Precision Test Worksheet 96-Tube Manifold Annual Buffer Switching Test	st211
Dispense Precision Test Worksheet 192-Tube Manifold	212
Evacuation Efficiency Test Worksheet 192-Tube Manifold	213
Dispense Precision Test Worksheet 192-Tube ManifoldAnnual Buffer Switching	
Test	
Evacuation Efficiency Test Worksheet 128-Tube Aspirate Manifold	
Peri-pump Dispense Precision and Accuracy Tests	
Dispense Precision and Accuracy Specifications	
Peri-pump Dispenser Test Materials	
Peri-pump Dispenser Precision and Accuracy Test Solutions	
Perform the Peri-Pump Precision and Accuracy Tests	
Syringe Dispenser Liquid Tests	
Syringe Dispenser Test Materials	
Syringe Dispenser Test Solutions	
Perform the Syringe Dispense Precision & Accuracy Test	221
Syringe Dispenser Precision & Accuracy Test WorksheetTest 1 / 96-Well Microplate/40 µL Dispense	223
Syringe Dispenser Precision & Accuracy Test WorksheetTest 1 / 384-Well Micro-	220
plate / 20 µL Dispense	224
Syringe Dispenser Precision & Accuracy Test WorksheetTest 2 / 96-Well Micro-	
plate/160 µL Dispense	225
Syringe Dispenser Precision & Accuracy Test WorksheetTest 2 / 384-Well Micro-	
plate / 80 µL Dispense	226
Syringe Dispenser Precision & Accuracy Test WorksheetTest 3 / 1536-Well Micro-	
plate / 6 µL Dispense	227
Syringe Dispenser Precision & Accuracy Test WorksheetTest 4 / 96-Well Micro-	220
plate/300 µL Dispense	
Troubleshooting	
Troubleshooting	
Error recovery	
Error Codes	
General Operation Troubleshooting	
Startup	
Microplate Scratches	
Cell-based Assavs	231

7

Washer Troubleshooting	
Fluid Aspiration	
Fluid Delivery	
Fluid Leakage	
Washer Manifold Movement	
Vacuum Pump Problems	
Peri-pump Troubleshooting	
Syringe Dispenser Troubleshooting	
Fluid Delivery	237
Fluid Leakage	
Dispense Manifold Movement	
Microplate Carrier Movement	
Syringe Movement	
Communication Errors	240
Appendix A: Peri-pump Test Worksheets	.242
Peri-pump Dispense Precision & Accuracy Test Worksheet1 µL Test for 1 µL Cas-	
sette	243
Peri-pump Dispense Precision & Accuracy Test Worksheet10 µL Test for 1 µL Cas-	
sette	244
Peri-pump Dispense Precision & Accuracy Test Worksheet1536 Test	245
Peri-pump Dispense Precision & Accuracy Test Worksheet5 µL Test	246
Peri-pump Dispense Precision & Accuracy Test Worksheet50 µL Test for 5 µL Cas-	
sette	
Peri-pump Dispense Precision & Accuracy Test Worksheet 10 μL Test for 10 μL Ca	
Peri-pump Dispense Precision & Accuracy Test Worksheet100 µL Test	249
Appendix B: Specifications	250
Physical Specifications	251
Performance Specifications	253
Washer	253
Peri-Pump Dispenser	253
Cassette Expected Lifetime	254
Syringe Dispensers	254
Appendix C: Error Codes and Conditions	256
Error Codes Overview	
When an Error Code or Error Message Appears:	
System Error Codes	
406 FX-I HC Software Error Codes	266

Appendix D: Changing the Washer Manifold	272
Changing the Wash Manifold	273
To reconfigure the instrument, first update the washer manifold setting:	274
Storing the Unused Manifolds	274
Appendix E: Safety Information	275
Safety Notices	276
Warnings and Precautions	
Electrical Hazards	278
Chemical/Environmental	280
Components	287
Intended Product Use	290
Chemical Compatibility Chart	294
In This Book	298

2 Introduction

This chapter introduces the 406 FX and provides instructions for contacting Technical Support.

Introducing the 406 FX Microplate Washer Dispenser	15
Package Contents	16
Optional Accessories	
Peri-pump Optional Accessories	
Features of the 406 FX	
BioStack Compatibility	24
406 FX Dispenser Comparison	

2 Introduction

Introducing the 406 FX Microplate Washer Dispenser

Introducing the 406 FX Microplate Washer Dispenser

The 406 FX offers up to three devices in one instrument: a microplate Washer, a peristaltic pump dispenser (the "Peri-pump"), and a dual Syringe pump dispenser. One or both dispensers are provided with every 406 FX.

One model is available: 406FX--SN.




Figure 2-1: 406 FX Washer Dispenser

	Device/Component	Description
1	Syringe Dispensers	Syringe-pump dispensers, each with an 8-, 16- or 32-tube manifold. Up to two syringe pumps may be connected.
2	Peri-pump Dispensers	Peristaltic, 8-channel dispensers with an entirely visible fluid path. Shown with optional secondary Peri-pump.
3	Dispense arm	Holds the Peri-pump Tip Holder and the Syringe dispenser manifolds.
4	Washer manifold	96- or 192-tube devices aspirate and dispense fluid; or 128-tube manifold aspirates fluid from 1536-well plates.
5	Touchscreen	On-board instrument control.
6	Priming trough	Waste reservoir for collecting priming fluid for the Peri-pump and Syringe dispensers.
7	Ultrasonic cleaning and priming trough Mist shield	All models support AutoClean ultrasonic cleaning of the washer.
8	Plate Carrier	Holds standard microplates for processing.

Package Contents

Description	Part number
Power supply, 24 V DC	02285
External Vacuum Pump Switch Box (with control cable)	2120505
USB cable	75108
Microplate carrier	2120508
Mist shield and 2 thumbscrews	2122074
Reagent bottle holder	7212034
Reagent bottle holder step, 2 count.	7212035
Phillips screwdriver for removing dispenser arm shipping bracket	01188
Syringe priming trough	2122043
Peri-pump priming trough, 2	7182043
406 FX User Manual, on USB flash drive	2121000N

Accessory Pack

Description	Part number			
Hex wrench, 9/64" for removing washer manifold shipping bracket	48434			
Hex wrench, 1/16" for removing magnets from syringe manifolds and for adjusting manifold stop screw	48713			
Phillips screwdriver 0 3 mm OD, for Peri-pump cassette calibration	98268			
Disposable syringe, 10 cc Luer, for Peri-pump cassette maintenance	49919			
Stylus, .012 OD no color, for cleaning 192-tube dispense manifold	7772029			
Stylus, .018 OD yellow, for cleaning washer and Syringe manifold dispense tubes, and 128-tube aspirate manifold tubes	7772030			
Stylus, .031 red, for cleaning washer manifold aspirate tubes	7772033			
Manifold stop height set tool for adjusting manifold stop screw	2122136			
Shipping set	2123008			

Optional Accessories

General Instrument Accessories

Description		PN
Agilent liquid testing solutions for instrument qualification tests	Wetting Agent Blue Test Dye	7773002 7773001
Liquid Handling Control Software 3.0		5330209

2 Introduction

Optional Accessories

Washer Accessories

Description	Part Number
Complete dispense/waste system 4 L bottles	2120531
Complete dispense/waste system High Flow 4 L bottles	2120532
Complete dispense/waste system High Flow 10 L bottles	2120533
Complete dispense/waste system High Flow 20 L bottles	2120534
Vacuum System Complete Assembly 115V/230V, 4 L bottles	7100746S
High Flow Vacuum System Complete Assembly 115V/230V, 4 L bottles	7100753S
Complete Direct Drain Waste System, 115V/230V, 4L bottles	1170560S
Complete Direct Drain High Flow Waste System 115V/230V, 4 L bottles	1170562
Direct drain waste system upgrade kit, 115V/230V, 4 L bottles (For use with any complete waste system)	1170561S
Complete dispense system. 4 L, Assembly	2120535
Complete dispense system.10 L, Assembly	2120536
Complete Waste System 115V/230V, 4 L bottles	7100746S
High Flow Vacuum System, Complete Assembly 115V/230V, 4 L bottles	7100753S
Vacuum Set with Level Sense Assembly, 4 L	7100543S
Vacuum Set with Level Sense Assembly, 10 L	7100582S
Vacuum Pump Assembly 115V/230V	7103035
High Flow Vacuum Pump Assembly 115V/230V	7100754

Repackaging and Shipping Instructions

Repackaging and Shipping Instructions

Please read the information provided below before preparing the instrument for shipment.

- Contact Technical Support before returning equipment for service.
- Decontamination prior to shipment is required by the U.S. Department of Transportation regulations.
- If the instrument has been exposed to potentially hazardous material, decontaminate it to minimize the risk to all who come in contact with the instrument during shipping, handling, and servicing. The Maintenance section contains decontamination instructions.
- Install the shipping hardware (see next section).
- The instrument's packaging design is subject to change. If the instructions in this document do not apply to the packaging materials you are using, contact Technical Support for guidance.
- Be sure to use packaging materials supplied by the manufacturer. Other forms of commercially available packaging are not recommended and can void the warranty.
- If the packaging materials have been damaged or lost, or if the same set has been used more than four times, order replacements.

CAUTION

Shipping Hardware. All shipping hardware must be removed before operating the instrument and reinstalled before repackaging the instrument for shipment.

Prepare the instrument for shipping:

LHC	T	ouc	h screen
1	Select Tools > Instrument Utilities > Washer		Select Instrument > Next > Next > Other. Select Park Manifold.
2	Select the shipping bracket link under Service Functions .		

- 3 Remove the Peri-pump's cassette and the Syringe dispensers' manifolds.
- 4 Wrap rubber bands around the Peri-pump pump heads
- 5 Remove the plate carrier and put it in the accessories box.
- 6 Uninstall and repack separately the Syringe pump and its accessories, if applicable.
- 7 Slide the transport rail into position next to the dispense arm. It will be secured in place by one of the shipping brackets.


Reinstall the shipping hardware by reversing the steps described to Remove the Shipping Hardware.

NOTE

The instrument's packaging design is subject to change over time. If the instructions in this section do not appear to apply to the packaging materials you are using, please contact Agilent Technical Support for guidance.


Repackaging and Shipping Instructions

1 Repack the instrument accessories box bottom tray, placing each item into its designated slot.



Repackaging and Shipping Instructions


- 2 Repack the instrument accessories box top tray, placing item into its designated slot.
- 3 Place this tray in the box and cover with the tray cover.



- 4 Close the flaps and tape shut.
- 5 Prepare the main box for the instrument and accessories


Repackaging and Shipping Instructions

- 6 Place the instrument box (A) on a flat, sturdy surface and insert the bottom (B).
- 7 Insert the box insert (C) so that the handholds of the insert align with the outer box. This insert will not fill the entire space.
- 8 Place the prepared accessories box between the insert and the box as see below.



Repackaging and Shipping Instructions

9 Place the bottom shipping foam insert into the box and then place the bottom into the foam insert. These support the instrument during shipping.






The instrument is heavy. Have another person help you with the next packing steps to avoid damage to the instrument.

Repackaging and Shipping Instructions

- 10 Place the instrument in the poly bag.
- 11 Carefully place the instrument in the box, fitting it into the bottom foam insert.
- 12 Place the top foam insert on top of the instrument and then place the cardboard over it.




13 Securely tape up the box for shipment.


Repacking the Syringe Dispenser

Repacking the Syringe Dispenser

After preparing the instrument for shipping, and reversing the installation steps, pack the Syringe dispenser as shown here:



Fill the foam tray with the Syringe dispenser accessories and place the tray on top of the inner foam box containing the pump and supply bottles.



5 Maintenance

Properly maintaining the 406 FX is the key to reliable performance.

Preventative Maintenance Overview	158
Recommended Maintenance Schedule	
Daily Maintenance	
Overnight/Multi-Day Maintenance	
Soak the manifold tubes in cleaning fluid	
Removing Protein Residuals and Fungi Growth	166
Periodic Maintenance	167
Washer Maintenance	170
AutoClean the Washer	170
Remove and clean the washer manifold	
Peri-pump Dispenser Maintenance	
Syringe Dispenser Maintenance	178
Decontamination	183
Long Shutdown (Prepare for Storage or Shipment)	188
Replace Components	189
Calibrate the Backlash for Syringe Dispenser	192

5 Maintenance

Removing Protein Residuals and Fungi Growth

Removing Protein Residuals and Fungi Growth

CAUTION

Bovine Serum Albumin. Solutions containing proteins, such as bovine serum albumin (BSA), will compromise instrument performance over time unless a strict maintenance protocol is adhered to.

When using protein solutions or similar fluids, Agilent recommends performing the following procedures to thoroughly flush out protein particles and other contaminants from the fluid path.

Also note some components can be autoclaved to sterilize them.

NOTE

Four-liter volumes specified in the following are approximate amounts.

Daily Practice with buffer and deionized water:

If the 406 FX is idle between plates for longer than 45 minutes, flush the proteins using buffer. Use buffer for the AutoPrime, as well, for immediate reuse of the instrument. Otherwise, if the tubes will be allowed to dry, flush the system with DI water to remove the buffer:

- 1 Fill a supply bottle with buffer solution. Connect the bottle to the Buffer A port or to the syringe dispenser.
- 2 Run the applicable DAY_RINSE protocol.
- 3 Enable AutoPrime for 60-minute intervals.

NOTE

Flush the buffer from the system using DI water when the device will be idle for an extended period, that is, repeat steps 1 and 2 using water.

At the end of the day:

- 1 Fill a supply bottle with buffer for each device (washer or syringe). Connect the bottle to the Buffer A port for the washer or to the appropriate Syringe dispenser port.
- 2 Run the applicable DAY_RINSE protocol three times.
- **3** Fill a supply bottle with deionized water and repeat steps 1 and 2.
- 4 Perform your regular Overnight/Multi-Day Maintenance routine.

Weekly or As Needed, use NaOH and HCl to remove proteins:

- 1 Flush the system with 0.1-0.5 N* NaOH (sodium hydroxide), followed by neutralization with an equivalent normality (0.1-0.5 N) of HCl (hydrochloride).
- 2 Rinse well with deionized water to remove the HCl.
- 3 Run the applicable W-DAY_RINSE protocol three times with deionized water if you plan to use the device immediately.

NOTE

* N = Normal solution, which contains 1 'gram equivalent weight' (gEW) of solute per liter of solution. The gram equivalent weight is equal to the molecular weight expressed as grams divided by the 'valency' of the solute.

Periodic Maintenance

Alternatively, use an Enzyme-Active Detergent:

- 1 Mix an enzyme-active detergent according to the manufacturer's directions to fill two four-liter supply bottles.
- 2 Run the W-DECONTAMINATE or S-DECONTAMINATE protocol, as appropriate, based on the syringe channel.

NOTE

Follow the on-screen prompts and messages. **Delay** messages contain instructions and information to help you walk through the procedure.

- 3 Respond to the **Delay** message and when ready, press **Continue**.
- 4 When the protocol is completed, determine if the appropriate Decontaminate protocol needs to be run again and repeat as necessary.
- 5 When the process is completed, connect a bottle containing four liters of deionized water and run W-DAY_RINSE or S-DAY_RINSE as appropriate, based on the syringe channel, three times to flush the system.
- **6** Repeat the procedure for the other Syringe dispensers as appropriate.

Direct Drain Only: Clean the Waste Bottle Filter

Periodically clean the direct-drain waste bottle filter (part number 01310):

WARNING

Direct Drain Waste Bottle Filter: Handle with care! When cleaning the waste bottle filter, be aware of any potential biohazards and handle appropriately.

- 1 Open the bottle and lift the cap and its tubing up and out of the bottle
- 2 Remove the fluid filter at the bottom of the tubing.
- **3** Wash the filter with hot water and a soft-bristle brush, if necessary.
- 4 Rinse the filter and reinstall it.



Periodic Maintenance

Periodic maintenance involves cleaning the components on a regular basis to keep the instrument running efficiently and in compliance with performance specifications. The recommended **frequency for cleaning components** is *at least monthly*. The risk and performance factors associated with your assays may require that some or all procedures be performed more frequently.

WARNING

Internal Voltage. Always turn off the power switch and unplug the power supply before cleaning

5 Maintenance

Periodic Maintenance

the outer surface of the instrument.

WARNING



Potential Biohazards. Some assays or specimens may pose a biohazard. Adequate safety precautions should be taken as outlined in the assay's package insert. Always wear safety glasses and appropriate protective equipment, such as chemical-resistant rubber gloves and apron.

WARNING

Liquids. Avoid spilling liquids on the instrument; fluid seepage into internal components creates a potential for shock hazard or instrument damage. If a spill occurs while a program is running, stop the program and turn off the instrument. Wipe up all spills immediately. Do not operate the instrument if internal components have been exposed to fluid.

CAUTION

Liquids. Do not immerse the instrument, spray it with liquid, or use a dripping-wet cloth on it. Do not allow water or other cleaning solution to run into the interior of the instrument. If this happens, contact Technical Support. Do not soak the touch screen.

CAUTION

Lubricants. Do not apply lubricants to moving parts. Lubricant on components in the carrier compartment will attract dust and other particles, which may cause the instrument to produce an error. Do not apply lubricants to manifold o-rings, channel-end seals, bottle cover seals, any tubing connection, or any surface that is a part of the fluid path.

CAUTION

Chemical Compatibility. Some chemicals may cause irreparable damage to the instrument. The following chemicals have been deemed safe for use in the instrument: buffer solutions (such as PBS), saline, surfactants, deionized water, 70% ethyl, isopropyl, or methyl alcohol, and 20% sodium hydroxide. Never use acetic acid, DMSO, or other organic solvents. These chemicals may cause severe damage to the instrument.

CAUTION

Caution: Bovine Serum Albumin. Solutions containing proteins, such as bovine serum albumin (BSA), will compromise instrument performance over time unless a strict maintenance protocol is adhered to.

Perform these maintenance tasks regularly:

- · Clean the Bottles below.
- Clean the Plate Carrier on the facing page.
- Peri-pump Dispenser Maintenance on page 175.
- Syringe Dispenser Maintenance on page 178.

Clean the Bottles

- Clean and rinse the supply bottles with deionized water before the first use, before each refill, and as necessary, to prevent bacteria growth.
- Empty the waste bottle often (at least daily), and firmly seat the waste bottle stopper.
- Rinse the covers every time the wash or rinse bottles are filled.
- Accumulated algae, fungi, or mold may require decontamination.
- Check the hex nuts securing the quick-disconnects to the bottle cap to ensure they are not loose or corroded.

NOTE

To ensure that fluid does not back up into the vacuum pump during operation, always operate the instrument with the **waste sensor cable** installed and the **waste detection sensor** enabled (the

Periodic Maintenance

sensor is enabled by default).

NOTE

If fluid collects in the overflow bottle, thoroughly rinse the level-switch assembly and bottle.

Clean the Plate Carrier

If liquid has overflowed onto the plate carrier, transport rail, or glide strips, some buildup may occur and prevent the microplate from seating correctly on the carrier. This buildup can interfere with plate transport. Weekly cleaning is recommended.

- 1 Turn the instrument off.
- 2 Lift the carrier up and off the transport rail.
- 3 Clean the carrier, rails, and glide strips using mild detergent and hot water, 70% isopropyl alcohol, or ethanol. Clean the priming trough, too.
- 4 If detergent was used, wipe the components with a cloth moistened with water. Use a clean, dry cloth to dry the components.
- **5** Reinstall the carrier:
 - Place it on the transport rail so the slot on its base fits into the **Y-axis Carrier Arm**.
 - If necessary, release the spring-loaded microplate clamp in the back left corner of the carrier to level the carrier on the base.

Clean the Exterior Surfaces and Mist shield

- 1 Turn off the instrument and disconnect the power cable.
- 2 Moisten a lint-free disposable towel with water or with water and mild detergent. **Do not soak the cloth**.
- 3 Remove the mist shield if it is attached. Wipe the inside and outside surfaces of the mist shield with the towel. Wipe the top surface of the instrument base and all exposed surfaces of the instrument.
- 4 If detergent was used, wipe all surfaces with a water-moistened cloth.
- 5 Use a clean, dry cloth to dry all wet surfaces.

Autoclavable Components

Autoclaving is an efficient method of sterilizing instrument components. It is a good alternative to some of the decontamination procedures for qualified items.

Do autoclave:	Do NOT autoclave:						
Peri-pump cassette*	Washer manifolds						
16-Tube and Low-density plate-type-specific Syringe dispenser manifolds	Plate carrier						
Priming trough inserts	32-Tube Syringe dispenser manifolds						
Syringe pump head (glass/stainless steel)	8-Tube Syringe dispenser manifolds (gray block holds two manifolds, 16 tubes total)						



Chemical Compatibility

For your own safety and to ensure the instrument's long life, follow these guidelines when choosing compounds to use.

Important: Avoid continuous contact with harsh chemicals. A best practice is to rinse all components with deionized/distilled water after contact with any strong acid, base, or solvent.

Materials listed as 1-20 in **Table 1** are cross-referenced as headings 1-20 in **Table 2**. Table 1 indicates which component contains the material using **P** for Peri-pump, **S** for the dual Syringe dispensers, and **W** for the washer module.

Table 1: Material Where Used

#	Material		Where Used
1	304 Stainless Steel	W	Inlet screen, feeder tube to manifold, vacuum switch, fluid pump, direct-drain inlet screen
2	316 Stainless Steel	Р	Steel dispense tips, pump rollers, carrier rails, cassette bolts
		S	Dispense tubes, syringe piston, cylinder head and fittings, check valve spring, aspirate tubes, aspirate arm release button
		W	Fluid inlet filter, feeder tube to manifold, vacuum switch, fluid pump, and vacuum filtration pressure regulator components
3	Aluminum		Instrument main structural plate
		Р	Center Holder, Tube Tensioner, pump head, dispense arm
4	Borosilicate Glass	S	Syringe cylinder, supply bottles (autoclavable models only)
5	Brass	Р	Threaded inserts in Center Holder and Tip Holder
6	EP (Ethylene Propylene)	S	Check-valve o-rings (organic solvent tubing set)
		W	Inlet valve, vacuum switch
7	ETFE (Ethylene	S	Manifold fitting
	tetrafluoroethylene)		
8	Neoprene	W	Manifold channel-end seals
9	Nylon	W	Inlet fitting, vacuum switch adjustment screw, carrier leveling feet
9	Polyethylene	SW	Supply Bottle
10	Polycarbonate		Plate carrier
		Р	Plate positioning dog, peristaltic pump cover latch
		W	Vacuum switch
11	Polyethylene	W	Buffer bottle
12	Polypropylene	Р	Molded dispense tips
		S	Bottle cap and fittings, check valves (organic solvent tubingset)
		W	Outlet fitting, fittings in bottles, inline fittings, float ball, vacuum filtration pressure gauge components, direct-drain module fittings
13	Polystyrene		Assay plates, mist shield
		Р	Peristaltic pump cover, RFID antenna cover
		W	Flow sensor
14	PTFE (polytetrafluoroethyler	ne) S	Syringe seal
		W	Optional check valves (PN: 68098) for fluid pump
15	CPVC (Polyvinyl chloride)	W	Inlet screen, inlet valve, manifold, waste sensor, ball retainer, vacuum filtration carrier, system orifice, and pressure regulator, direct drain bottle sensor



#	Material		Where Used
16	Ryton PPS	Р	Tube Organizer, priming trough, priming trough insert, RAD 8-to-1 chute
		S	Manifold, manifold plugs, priming trough insert
		W	Outlet valve, fluid pump, inlet valve, direct-drain fluid pump
17	Santoprene	W	Fluid pump, direct-drain pump check valves
18	Silicone	Р	Tubing
		S	Tubing, manifold o-rings, check valve o-rings
		W	Inlet tubing, outlet tubing, o-rings, vacuum filtration pressure regulator tubing and o-rings, direct-drain module tubing
19	Ultem (polyetherimide)	Р	Center Holder, Tube Tensioner, Tip Holder
		S	Check-valves
		W	Outlet valve, cell wash inlet valve, vacuum filtration valve, direct-drain valve
20	Viton	W	Outlet valve, direct-drain valve



Table 2: Chemical Compatibility

Chemical		1-20 represent Materials listed above in the Where Used List																		
	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
Key					ne)	(a)														
A - No effect	<u>a</u>	_	lass		pyle	ene /lene	<u></u>	ate	ЭС	au e	ā	nted Ilene	ated ride)				ated ride)		ide)	
B - Slight effect	Stee	Aluminum	Borosilicate Glass	Brass	(Ethylene Propylene)	ETFE (Ethylene tetrafluoroethylene)	Noryl (PPO)	Polycarbonate	Polyethylene	Polypropylene	Polystyrene	FEP (fluorinated ethylene propylene)	CPVC (chlorinated polyvinyl chloride)	Ryton	Silicone	Ultem	CPVC (chlorinated polyvinyl chloride)	Ryton	Ultem (polyetherimide)	Viton
C - Moderate effect	316 S.	lum	Silics	Bra	ylene	E (E	oryl	ycar	lyetl	ypro	olyst	(flu	C (ch vinyl	Ry	Silic) j	C (ch vinyl	Ryt	Ulte	Vit
D - Severe effect	ω,	٩	Boro		(Eth)	ETF	Z	Pol	Pc	Pol	ď	FEP	CPV(poly				CPV(poly		lod)	
N - No data					EP	te						O								
Acetic Acid, 5%	Α	В	Ν	D	Ν	А	А	А	А	А	D	А	D	А	А	А	D	Α	А	А
Acetic Anhydride	А	А	N	D	А	N	D	D	С	В	D	А	D	Α	А	N	D	А	А	N
Acetone*	D	D	D	D	D	D	D	D	D	D	D	D	D	D	D	D	D	D	D	D
Acetonitrile	Α	В	А	Ν	С	А	N	D	A	А	D	А	D	Α	D [†]	D	D	Α	D [†]	D
Ammonia 10%	А	А	N	Ν	Ν	N	А	D	N	A	В	A	В	Α	D	D	В	А	D	D
Benzyl Alcohol	В	В	Α	Ν	В	А	D	N	А	А	D	А	D	Α	А	N	D	Α	А	N
Chloroform	А	В	А	В	D	A	D	D	D	С	D	А	D	Α	D	D	D	Α	D	D
Detergents 1%	А	В	А	Ν	А	А	А	А	А	A	А	А	А	Α	А	А	А	Α	А	А
Dimethylformamide	В	А	Α	Ν	В	А	D	D	С	Α	D	А	D	Α	А	N	D	А	А	N
DMSO (Dimethyl sulfoxide)	А	А	А	N	В	С	В	D	А	А	D	А	D	Α	C [†]	D	D	Α	C [†]	D
Ethyl Alcohol 70%	А	А	Α	A	А	А	А	В	В	А	А	А	В	Α	В	А	В	А	В	А
Ethylene Oxide	В	D	A	D	С	А	А	С	А	D	С	А	С	D	А	N	С	D	А	N
Formaldehyde 37%	Α	В	А	А	Α	А	А	А	D	А	N	А	Α	Α	С	А	Α	А	С	А
Hexane	Α	А	А	А	D	A	В	D	А	В	D	А	В	Α	D	А	В	А	D	А
Hydrochloric Acid 20%	D	D	Α	D	А	A	А	В	А	В	С	А	А	D	D	А	Α	D	D	А
Hydrogen Peroxide 10%	В	А	Α	D	А	А	А	А	А	А	А	А	Α	С	А	А	Α	С	А	А
Isopropyl Alcohol 70%	Α	А	Α	Α	А	А	А	А	В	А	А	А	В	Α	А	А	В	А	А	А
Methyl Alcohol 70%	А	А	Α	А	А	А	А	В	А	А	N	А	Α	Α	D	А	Α	Α	D	А
Methylene Chloride	В	С	Α	Α	Ν	А	D	D	D	В	D	А	D	Α	D	D	D	А	D	D
Phosphoric Acid 10%	D	С	N	D	А	А	А	А	А	А	В	А	В	Α	С	А	В	Α	С	А
Propylene Glycol	В	В	N	Ν	Ν	N	N	В	А	А	А	А	С	Α	А	N	С	Α	А	N
Sodium Chlorate	В	С	N	Ν	Z	А	А	А	N	А	N	А	А	А	С	N	А	А	С	N
Sodium Hydroxide 20%	В	D	А	В	В	А	А	А	А	А	А	А	А	Α	А	А	А	А	А	А



Chemical	1-20 represent Materials listed above in the Where Used List																			
	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
Sodium Hypochlorite <20%	С	D	Α	D	В	А	А	С	А	А	А	А	А	C	В	В	А	С	В	В
Sodium Hypochlorite 0.50%	В	D	А	D	В	А	А	С	Α	А	А	А	A	С	В	А	А	С	В	А
Sulfuric Acid <10%	В	В	А	N	А	С	А	А	Α	А	А	А	А	А	С	А	А	А	С	А
Trichloroethylene	В	D	N	N	N	А	D	N	D	С	D	A	D	А	D	D	D	А	D	D
Virkon 10%	А	D	А	N	А	N	А	А	А	А	А	N	А	А	А	N	А	А	А	N

^{*} Acetone is not recommended for use with Agilent instruments.

† Exposure to DMSO and Acetonitrile may cause the silicone tubing to swell, increasing the volume of fluid Dispensed. The magnitude of this effect will vary with concentration and exposure. Re-calibration of the cassette may be required.